1 mL Luer Lock Syringe REF QIT001
3 mL Luer Lock Syringe REF QIT002
20 mL Luer Lock Syringe REF QIT003

STERILE R

CE 2019
0459

Guerbet |

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LUER LOCK SYRINGE - PRECAUTIONS FOR USE

EN

Single use syringe, intended for injecting fluids into the body, manually, during procedures using Lipiodol® Ultra Fluid

• This device must be used by qualified healthcare professionals. • This device does not contain latex. • Single-use device: this device is to be used once on one patient only. • If reused, single-use devices may no longer function effectively or perform as intended over time. The patient may also be at risk of contamination if the device is reused. • Before use, always check that the medical device and its packaging are undamaged. • Before and during use, ensure the safety, hygiene and asepsis instructions are followed. • For safe & proper use of the product, please carefully read & follow the instruction(s) for use of the drugs used with the device. Ensure that the use of this kit respects the valid rules of good practice & the marketing authorization (MA) of the used products. • This device is for manual use. • Connections are compatible with Luer 6%. • All Luer connections should be gently hand tightened without over tightening to ensure their security and to prevent damage to the device. • Apply axial controlled pressure on the syringe's plunger to avoid uncontrolled injection and damage to the device. • After use, the devices must be treated as hazardous health care waste, according to the regulations in force for the device and the administered fluids. • If the used medical device is to be returned to the manufacturer, ensure it is decontaminated first, according to the regulations in force for the device and the administered fluids.

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
4-WAY STOPCOCK - PRECAUTIONS FOR USE


EN

Single use device intended for fluid flow directional control and for providing access ports for administration of Lipiodol® Ultra Fluid

• This device must be used by qualified healthcare professionals. • This device does not contain latex. • Single-use device: this device is to be used once on one patient only. • If reused, single-use devices may no longer function effectively or perform as intended over time. The patient may also be at risk of contamination if the device is reused. • Before use, always check that the medical device and its packaging are undamaged. • Before and during use, ensure the safety, hygiene and asepsis instructions are followed. • For safe & proper use of the product, please carefully read & follow the instruction(s) for use of the drugs used with the device. Ensure that the use of this kit respects the valid rules of good practice & the marketing authorization (MA) of the used products. • This device is for manual use. • Connections are compatible with luer 6%. • All Luer connections should be gently hand tightened without over tightening to ensure their security and to prevent damage to the device. • After use, the devices must be treated as hazardous health care waste, according to the regulations in force for the device and the administered fluids. • If the used medical device is to be returned to the manufacturer, ensure it is decontaminated first, according to the regulations in force for the device and the administered fluids.

!! This device is not to be used with surgical glue.

!! To disconnect the stopcock, hold the downstream device then unscrew the mobile ring without turning the stopcock.

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